Opto-Electronics

Features

- Single phase instrument
- · Differential or single-ended input
- Gain settings from 3 µ V to 1V
- High performance wide bandwidth input gain stage
- Analogue meter for display of output signal
- Output Offset controls
- Output time constants from 100 μs to 30s
- 1F and 2F reference signal operation
- 90° step and fine phase control



Input Signal Channel

Input	High or low impedance differential or single ended via front panel BNC
Sensitivity	3 μ V to 1V (for 1V output) switched in 1, 3, 10 steps
Input Impedance	10 ¹² W//1nF, dc coupled
Frequency	10Hz to 100kHz
Maximum Inputs	± 16 V before input protection circuitry comes into operation.
	The input BNC has been tested for electrostatic discharge damage.
Gain Accuracy	1%
Gain Stability	200ppm/℃
Dynamic Reserve	60dB limited by a maximum signal input noise voltage of 10V

Demodulator

The output of the signal input stage is processed using a very high bandwidth demodulator to recover the input signal. Offsets introduced at this stage are automatically removed via novel feedback mechanisms.

Low Pass Filter

Time Constant	100 µs to 30s in 1, 3, 10 steps
Output	$\pm 1 V$ output corresponds to full scale Input. Short circuit protection included.
Offset	Up to 1x full scale, switchable on or off

Reference Channel

Frequency		10Hz to 100kHz
Input Impedance		$5.6 M_{\Omega}$ ac coupled
Trigger	-Sine	100mV rms min (15V max.)
	-Pulse	5V, 95% mark/space ratio min.
Acquisition time		10s max.
Phase control		90° steps + fine shift in range 0° - 100°
Phase Drift		0.1° /℃

General

Power	115Vac, 230Vac; 50-60Hz; 10VA max.
Mechanical	440mm (W) x 87mm (H) x 190mm (D)
Temperature range	0-50℃ (operational)

Ordering Information

Model	Descriptions
410	A complete cased analogue single phase lock-in amplifier with IEC power lead.
410-RMK	A complete cased analogue single phase lock-in amplifier with IEC power lead
	but with lid designed for 19" rack mounting.